

o858:medium precipitation, granular, retention range approx. 7-12 µm
medium fast

- Detection of metals in milk fat
- General filtration paper, germination test
- For the reliable clarification of viscous liquids
- Creped filter papers have a particularly large surface area and correspondingly shorter filtration times



ITEM NO.	DESCRIPTION	DIA-METER	SORT	QUAN-TITY	UNIT
HMDP5891047	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	47 mm	589/1	100	piece
HMDP5891055	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	55 mm	589/1	100	piece
HMDP5891070	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	70 mm	589/1	100	piece
HMDP5891090	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	90 mm	589/1	100	piece
HMDP5891110	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	110 mm	589/1	100	piece
HMDP5891125	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	125 mm	589/1	100	piece
HMDP5891150	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	150 mm	589/1	100	piece
HMDP5891240	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	240 mm	589/1	100	piece
HMDP58920125	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	12,5 mm	589/2	1000	piece
HMDP58920127	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	12,5 mm	589/2	1000	piece
HMDP5892047	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	47 mm	589/2	100	piece
HMDP5892055	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	55 mm	589/2	100	piece
HMDP5892070	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	70 mm	589/2	100	piece
HMDP5892090	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	90 mm	589/2	100	piece
HMDP5892110	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	110 mm	589/2	100	piece
HMDP5892125	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	125 mm	589/2	100	piece
HMDP5892150	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	150 mm	589/2	100	piece
HMDP5893045	Type 589/3, quantitative filter paper, slow, 84 g/qm, filter circles	47 mm	589/3	100	piece