

**o858:**medium precipitation, granular, retention range approx. 7-12 µm  
medium fast

- Detection of metals in milk fat
- General filtration paper, germination test
- For the reliable clarification of viscous liquids
- Creped filter papers have a particularly large surface area and correspondingly shorter filtration times



ITEM NO.	DESCRIPTION	DIA-METER	SORT	QUAN-TITY	UNIT
<b>HMDP5891047</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	47 mm	589/1	100	piece
<b>HMDP5891055</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	55 mm	589/1	100	piece
<b>HMDP5891070</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	70 mm	589/1	100	piece
<b>HMDP5891090</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	90 mm	589/1	100	piece
<b>HMDP5891110</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	110 mm	589/1	100	piece
<b>HMDP5891125</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	125 mm	589/1	100	piece
<b>HMDP5891150</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	150 mm	589/1	100	piece
<b>HMDP5891240</b>	Type 589/1, quantitative filter paper, fast, 79g/qm, filter circles	240 mm	589/1	100	piece
<b>HMDP58920125</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	12,5 mm	589/2	1000	piece
<b>HMDP58920127</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	12,5 mm	589/2	1000	piece
<b>HMDP5892047</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	47 mm	589/2	100	piece
<b>HMDP5892055</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	55 mm	589/2	100	piece
<b>HMDP5892070</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	70 mm	589/2	100	piece
<b>HMDP5892090</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	90 mm	589/2	100	piece
<b>HMDP5892110</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	110 mm	589/2	100	piece
<b>HMDP5892125</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	125 mm	589/2	100	piece
<b>HMDP5892150</b>	Type 589/2, quantitative filter paper, medium fast, 85g/qm, filter circles	150 mm	589/2	100	piece
<b>HMDP5893045</b>	Type 589/3, quantitative filter paper, slow, 84 g/qm, filter circles	47 mm	589/3	100	piece